

IN-CUT



VECTEC

SINGLE USE
ENDOSCOPIC DEVICES





IN CUT

HIGH QUALITY SINGLE USE LAPAROSCOPIC SCISSORS



REF: CAB 3000
« Single use monopolar cable »



Reference	Working length	Diameter	Blade cut length	Connector
IC5 330	330 mm	5 mm	17 mm	Radial
IK5 330	330 mm	5 mm	17 mm	Axial
T5 330	330 mm	5 mm	12 mm	Radial
C5 470 *	470 mm	5 mm	17 mm	Radial
C5 250 *	250 mm	5 mm	17 mm	Radial

* box of 5 individually packaged units

Ergonomic handle allowing a precise surgical gesture.

IK5 330 model developed for gastrointestinal surgery.

IC5 330 model dedicated to gynecological surgery.

Double action curved blades facilitating the dissection.

Monopolar coagulation with extra long electrical insulation up to the active portion of the blades.

Packaging: box of 25 individually packaged units.



- Certified Management System
- EN ISO 9001
- EN ISO 13485
- ISO 13485 under CMDCAS



MD Class IIb
Made in France
See instructions for use

Graphisme : PoissonsLune - Novembre 2014
Non contractuel pictures.



VECTEC

Bioparc - Vichy - Hauterive
03270 Hauterive
FRANCE

Phone : +33 (0)4 70 58 41 70
Fax : +33 (0)4 70 58 41 79
Email : vectec@vectec.fr
http://www.vectec.fr